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EN | Instructions for use
One-piece, latex-free tube holder for securing endotracheal tubes (oral/nasal), double lumen tubes, laryngeal tubes and laryngeal masks. The additional adhesive hook and loop tape strengthens the tube holder. It is applied in adolescents and adults.

1. Useful life
Replace the product in case of contamination / defect / according to hygienic guidelines. Maximum operating life 24 hours.

2. Notice on reporting requirements
Please inform the manufacturer and the competent authority of the Member State in the event of any serious incident that has occurred in relation to the device.

3. Notice on medical devices
Medical device of Class I. Accessory product for endotracheal tubes, double lumen tubes, laryngeal tubes and laryngeal masks.

4. Contraindications
Do not use in case of injuries or accesses in the operating area of the tube holder or allergies to the materials used..




5. Warning notices
• The tube holder is applied on the connected tube or on the connected mask.

• Apply the tube holder in twos and constantly ensure a secure fit of the tube or mask.

• The safe hold of the fixation has to be checked regularly.

• If the fixation is too loose, extubation can be caused. A slight tension on the tube holder is optimal.

• Make sure that the patient's skin is not exposed to the hook and loop strip of the tube holder, as this may cause irritations and open wounds.

REF	Article
20-02	TuBo-Clip S
20-06	TuBo-Clip S, XL
20-99	TuBo-Clip S, Klebeklettstreifen

Patientenzielgruppe: Intubierter Patienten
Patient target group: Intubated Patients

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