Pylant Monitor | Pylant Monitor
Cuffdruckmesser | Cuff Pressure Manometer


Novo Klinik-Service GmbH
Technologiepark West - Haus 15
Zum Frenser Feld 1 | 50127 Bergheim | Deutschland | Germany
Tel. +49 2271 99404-0 | Fax +49 2271 99404-100
info@its-novo.de | its-novo.de

Inoris Medical Sp. z o.o. | Poland
58-530 Kowary | Poland | Tel. +48 663 977 566 | www.inoris.pl

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EN | INSTRUCTIONS FOR USE

GENERAL USE

The Pylant Monitor is a disposable endotracheal tube (ETT) cuff pressure manometer. It is used to indicate the pressure within the ETT cuff to prevent inadvertent over or under inflation of the endotracheal tube cuff. The device is used in adolescents and adults. Only medical professionals may use the product after prior product training.

NOTICE ON MEDICAL DEVICES

Medical device of Class I.

CONTRAINDICATIONS

N/A

USEFUL LIFE

Once the desired ETT cuff pressure is achieved the Pylant Monitor assembly can be removed or just the syringe can be removed leaving the ETT cuff at the desired pressure. The Pylant Monitor is air tight when connected to the inflation port of the ETT and may be left inline if determined safe by the practitioner.

SINGLE USE NOTICE

The product is a disposable product. The reuse of disposables can result in contamination and/or effects on the functioning.

NOTICES ON USE

Inspect the device and check for leaks before use.
Pylant Monitor is connected to the ETT or tracheal cannula inflation port (pilot balloon) at one end and a syringe at the other end.
For endotracheal tubes and tracheal cannulas, preferably use a 20 mL syringe.
For endotracheal tubes and tracheal tubes, the pressure should be between 20 cm – 30 cm H₂O; the arrow is in the green area of the scale.
If the blue pointer is in the yellow range, block the cuff using the syringe already connected until the optimum pressure is built up in the cuff and the pointer is in the green range. If the pointer moves into the red area, the cuff must be unblocked.

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