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Mobilitätstrainer

Art.-Nr. 11400 050

PZN 09338545

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Mobility trainer

Item No. 11400 050

Please read these notes carefully before use. Be sure to follow the instructions and keep them in a safe place! If this aid is to be used by other persons, they must be fully informed of these safety warnings before use.

This mobility trainer is meant to help the user to carry out regular mobility training for arms and legs. The number of rounds achieved and the length of the training time can be read from the display. The resistance can be raised or reduced with the aid of a screw, according to the degree of increase.

The training can thus be individually increased according to the user's own needs. The aim is to maintain and increase existing muscle power. To do this, it is recommended that 2-3 exercise phases of 5-10 minutes per day are carried out.

Safety Instructions

- Please pay attention to loose, damaged, worn out or missing parts during assembling and before usage!
- Please replace damaged parts immediately, because otherwise the function and safety is affected!
- Regularly check the tight fit of the screws and nuts!
- This training device is only suitable for use at home, not for professional use!
- You should consult your doctor before starting the exercise and in case of occurrence of indisposition during the work out!
- Children should only use the mobility trainer after being given exact instructions and under supervision!
- Comfortable clothing should be worn for doing the exercises. However, you should avoid wide trouser legs or sleeves, which may get caught by the pedals!
- Before you begin exercising ensure that the mobility trainer is standing on a firm, level and non-slip surface!
- Before you begin training it is recommended that you do several stretching exercises. Begin pedalling slowly and increase the speed until you have reached the desired performance level! Do not stop the exercises suddenly, but slowly reduce the movements until your pulse has reached its normal level once again!
- Never put your full bodyweight on the appliance! The appliance is not designed to support this!
- The metal parts may become hot during use. Allow these parts to cool after use and do not pick up the unit immediately by the pedal arm.

Scope of delivery

1 Base unit, 1 battery (LR 44, already inserted in the cycle counter), 1 resistance screw, 2 velcrostraps

Assembly Instructions

Remove all packaging material. Fold both legs of the mobility trainer outwards until the locking pins engage. Place it in front of you and turn the resistance screw in the thread on the opposite side of the display. The more you tighten the screw, the greater the degree of difficulty.

The loops on the pedals are already pre-mounted – they can be adjusted in two different widths. To do this, loosen the loops on the inside of the pedals and fasten them again to the inner or outer holes.

Display

It is inserted into the frame of the body and can be easily taken out for battery control or change. For that purpose you stick a screwdriver into the little gap at the side of the round counter and carefully lift it up. If necessary you can loosen the cable connector.

When reinstalling make sure that the cable connector is in the right place and that the round counter is correctly inserted.

By pressing the red start button or simply by moving the pedal, the round counter is activated. It offers the following display functions, which can be selected by repeatedly pressing the red button:

Scan: In this mode the display changes automatically between the different display functions, here you have got a successive overview of all work out results.

Time: In this mode the duration of the current work out is displayed. The maximal duration is 99:59 min.

Count:	In this mode you are shown the completed rounds of the current work out.
Calories:	In this mode the burnt calories of the current work out are displayed. This is only an approximate value.
Laps per minute/	This mode displays the speed in laps per minute.
R.P.M.:	Please replace the battery in case the display can not be read clearly anymore. After a disuse of approx. 4 minutes the round counter switches off automatically. To reset all data press red button for about 4 seconds.

Instructions for use

Set up the mobility trainer in front of you. Begin your training with a light resistance, loosening the screw to do this.

For the leg training, sit down on a stable object e.g. a chair or the side of the bed, and put the device in front of you on the ground. If the user cannot sit on the side of the bed, the mobility trainer can also be set up on the foot of the bed using a support. Make sure that the support cannot slip. You can fasten the mobility trainer to the bed-frame with the help of the velcro strips that are included.

For the arm training you should place the mobility trainer on a firm standing surface, e.g. a table.

To prepare the muscles for the work out you should do some warm and loosening-up exercises before you start.

To give your body the possibility of regeneration, you should include time for a short relaxation after the work out. In doing so you can e.g. keep on pedalling in the lowest level possible.

Cleaning

The body and the pedals can be cleaned with a damp cloth. Do not use any aggressive cleaning agents. Afterwards use a soft cloth for wiping. Make sure that no water enters the equipment. Do not allow moisture to penetrate the lap counter.

Storage

During disuse you should store the mobility trainer in a dry place. You can fold up the legs of the mobility trainer in order to store it easily. To do this, screw tight the pins with the key-rings and fold both legs. ☀️ ☔

Disposal

Batteries may not be disposed of in household waste, they have a crossed out waste bin marked on them – see illustration. Every consumer is legally obliged to give batteries to a collection point within the local area, or to a shop, so that they can be disposed in an environmentally friendly way.

Batteries that contain toxic elements are also marked with the chemical signs of the pollutants (Cd = cadmium; HG = mercury; Pb = lead).

Please consult your local disposal company for disposal of the movement trainer and lap counter.

Technical data

- Length: approx. 47 cm
- Height: approx. 24 cm
- Width: approx. 37 cm
- Weight: approx. 2,5 kg

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